

SMARTER PRACTICO

The versatile operating table Practico has been designed to improve ergonomics with the industry's widest range of posture possibilities. It is powered by smart controls that help you work more fluently.



SUPER HIGH AND SUPER LOW

Improved ergonomics. Three different models, with motorized



STREAMLINED DESIGN

Human-centric design at its best. Perfect C-arm access. Easy to

leg sections as option.
The new smarter Practico is one of the lowest mobile tables in the market!

HIGH SWL

Stronger than ever. The new smarter Practico has improved safe working load and width extenders to help you to treat every patient regardless of their size and shape.

keep clean even in the busiest patient flow situations. No mechanical pedals.

SMART FEATURES

Practico is easy and intuitive to use. The user interface is similar to other Merivaara's products. Smart features also include for example safety sensors and memory positions.

FENNIAPRIZE 20

HONOURABLE MENTION

FENNIA PRIZE 2020 HONOURABLE MENTION

Smarter Practico® by Merivaara has been awarded the Fennia Prize Honourable Mention for the year 2020. The Fennia Prize is Finland's most prestigious business design competition.

Producto	Mesa de cirugía Practico		Mesa de cirugía Practico		Mesa de cirugía Practico	
Anchura (tope de la mesa)	540 mm	21,26 in	540 mm	21,26 in	540 mm	21,26 in
Ajuste de la altura	585 - 885 mm	23 - 34,8 in	640 - 1130 mm	25,2 - 44,5 in	540 - 1030 mm	21,3 - 40,6 in
Carga de trabajo segura (SWL)	280 kg	617,3 lb	280 kg	617,3 lb	280 kg	617,30 lb
Trendelenburg	30 °		30 °		30 °	
Trendelenburg inverso	30 °		30 °		30 °	
Ajuste del respaldo	-40 - 70 °		-40 - 70 °		-40 - 70 °	
Somier (secciones)	4-6		4-6		4-6	
Desplazamiento longitudinal de la parte superior de la mesa	390 mm	15,4 in	390 mm	15,4 in	390 mm	15,4 in
Inclinación lateral	-20 - 20 °		-20 - 20 °		-20 - 20 °	
Capacidad máxima de elevación	460 kg	1014,13 lb	460 kg	1014,13 lb	460 kg	1014,13 lb
Product Code	145000		145000		145000	