

# MAYO MESA DE INSTRUMENTOS

## Hydraulic instrument table

Hydraulic height adjustment 950–1335 mm

Lockable in all positions

Different tray top sizes

3 swivel castors 75 mm:n, 2 with brakes

*Tray size 600 mm x 400 mm*

Order number: I30021

*Tray size 700 mm x 500 mm*

Order number: I30024

*Tray size 850 mm x 600 mm*

Order number: I30027



Longitud min-max	600 - 850 mm
Ajuste de la altura	950 - 1335 mm
Tamaño de las ruedas	75
Product Code	I30021, I30024, I30027